

Art.-Nr.: KG1362
Ch.-B.: 023BM
Inhalt: 3 x 1 ml
Verw. bis: 2011-03

INTENDED USE

This product is intended for in-vitro use in the calibration of BGT Beta 2 microglobulin.

CHARACTERISTICS

Human source material from which this product has been derived has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE.

FDA approved methods have been used to conduct these tests.

However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

STORAGE & STABILITY

Unopened calibrator is stable until the expiry date printed on the product label when stored between 2°C and 8°C.

Once reconstituted the components are stable for 1 month at 2°C to 8°C or 3 months at -20°C.

PREPARATION FOR USE

Open the vial carefully, avoiding any loss of material and reconstitute with **1 ml** of distilled water. Replace the rubber stopper, mix gently and bring to room temperature for about 10mins before use. Ensure that all traces of dry material are dissolved.

MATERIALS REQUIRED BUT NOT PROVIDED

Distilled water

Volumetric pipette capable of accurate measurement of 1 ml

The assigned values are the means of at least 10 replicate determinations using clinical chemistry systems in a single laboratory, to the First @2-m International Reference from WHO.

Analyte	Target	
	mg/l	g/l
Beta 2 microglobulin	9	0.009

04 Dec '08