

Art.-Nr.: KG2492
Ch.-B.: 282IE – 287IE
Inhalt: 6 x 1 ml
Verw. bis: 2009-02

CHARACTERISTICS

The calibration material is prepared from purified Human IgE in a stabilised matrix.

VALUE ASSIGNMENT

Values have been assigned using an Immunoturbidimetric method standardized to the 2nd International Reference preparation 75/502.

SAMPLE PREPARATION

Calibrators are lyophilised. Reconstitute with 1ml of double distilled water and gently mix for 30 minutes, avoiding bubble formation.

STABILITY

The Calibrators are stable up to the expiry date when stored at +2 to +8°C.
After reconstitution the material is stable for 4 weeks at +2 to +8°C.

LOT NO.	IgE IU/ml	EXPIRY DATE
282IE	0.0	2009-02
283IE	52.5	2009-02
284IE	89.2	2009-02
285IE	201.1	2009-02
286IE	533.7	2009-02
287IE	1078	2009-02

For in vitro diagnostic use only, do not pipette by mouth, exercise the normal precautions required for handling laboratory reagents.

These products contain human material and should be handled as though capable of transmitting infectious disease.

This material contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.

Dispose of this material according to local regulations.

Revised 04 May '07