

**Chargen-B.: 863EC**  
**Inhalt: 12 x 5ml / 12 x 5ml**

**Artikel-Nr.: KG3107 / KG3108**  
**Verwendbar bis: 2011-11**

#### **INTENDED USE**

This product is intended for in vitro diagnostic use in the quality control of Immunoassays on clinical chemistry systems.

#### **SAFETY PRECAUTIONS AND WARNINGS**

For in-vitro diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

Human source material from which this product has been derived has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests.

However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

Health and Safety Data Sheets are available on request.

#### **STORAGE AND STABILITY**

**OPENED:** Store refrigerated (+2 to 8°C). Thawed serum is stable for 7 days at +2 to +8°C, with the following exceptions: Folate is stable for 3 days at +2 to +8°C. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.

**UNOPENED:** Store frozen at -20°C to -80°C. Stable to expiration date printed on individual vials.

#### **PREPARATION**

1. Allow the frozen control to thaw at room temperature (+18°C to 25°C) until completely thawed. Swirl the contents to ensure homogeneity.
2. Refer to the control section of the individual analyser application.
3. Refrigerate any unused material. Prior to reuse, mix contents thoroughly.

#### **MATERIALS PROVIDED**

Liquid Immunoassay Premium Level 3: 12 x 5ml Tri-Level: 12 x 5ml

#### **MATERIAL REQUIRED BUT NOT PROVIDED**

Volumetric Pipette

#### **ANALYTES**

AFP; Aldosterone; Amikacin;  $\epsilon$ 2 Microglobulin; Carbamazepine; hCG; Cortisol; DHEA-Sulphate; Digoxin; Estriol; Ethosuximide; Ferritin; Folate; Free T3; Free Thyroxine (FT4); FSH; Gentamicin; Growth Hormone; IgE; Insulin; LH; Oestradiol; Paracetamol; Phenobarbital; Phenytoin; Primidone; Progesterone; Prolactin; PSA (Free); PSA (Free); 17-OH-Progesterone; Salicylate; SHBG; Testosterone; Testosterone (Free); Theophylline; Total T3; Total T4; Tobramycin; TSH; T Uptake; Valproic Acid; Vancomycin; Vitamin B12.

**TYPICAL VALUES FOR THIS UNASSAYED CONTROL CAN BE OBTAINED AT:**  
**[www.biogentechnologies.de](http://www.biogentechnologies.de)**